

Selecting the appropriate VAP Shower Protection Pouch



The VAP Shower Protection Pouch comes in three sizes to accommodate your exit site dressing and your personal preferences.

To select the VAP pouch that is most suitable for you please use the template below to determine which VAP shower protection pouch you require. Your dressing should sit inside the aperture of the pouch, allowing a margin of approximately 1cm between the pouch and the dressing to ensure that the dressing is not inadvertently dislodged when your pouch is removed after showering.

Please seek advice from your Dialysis Unit if you are unsure which size you require.

All our Shower Protection Pouches are available on NHS prescription and may be obtained through the Independence Direct Home Delivery Service.

When you have selected the size you require please simply choose how you would like to order from the options below.



Call us on our Freephone number 0800 160 1399



Visit our website www.ipl.uk.com/irect and complete our online order form



Email us at enquiries@independencedirect.co.uk



Please ensure that you select "Actual Size" when printing this document to ensure that the template is accurate



ORDER FROM INDEPENDENCE DIRECT – 0800 160 1399
enquiries@independencedirect.co.uk – www.independencedirect.co.uk